

Updated 8/10

The Freedom[®] Knee components represent a true advancement in size and fit compatibility. The component sizing chart showcases the versatility of the system and the flexibility it provides surgeons in selecting the most appropriate option for each patient.

		Femoral Component Size									
		A	B	C	D	E	F	G	H		
Tibial Component Size (all-poly)	A1	■									
	B2		■								
	C2			■							
	C3			■							
	D2				■						
	D3				■						
	D4				■						
	E4					■					
	E5					■					
	F5						■				
G7							■				
Tibial Tray Size	1	■	■	■	■						
	2	■	■	■	■						1/2
	3			■	■	■	■				3/4
	4			■	■	■	■				5/6
	5					■	■	■			7/8
	6					■	■	■			
	7							■	■		
	8								■	■	
										Tibial Liner Size / Match	

NOTE




- The Freedom Total Knee[®] System’s Size “D” femoral component acts as a transitional size between the smaller “C” and larger “E” femoral sizes. Whereas the Freedom Knee’s components typically increase in the M/L dimension in 4mm increments, the “D” is only 2mm larger than the “C” and 2mm smaller than the “E” in the M/L dimension. This transitional size allows for more flexibility in sizing up or down between tibial and femoral components and provides a greater range of fit options in pairing these components.
- The above chart represents the core product sizes available for the Freedom[®] Knee. Additional sizes are available upon request.

COMPONENT	AP	ML
Femoral A (Universal)	51	54
Femoral B (Universal)	54	58
Femoral C (Universal)	58	62
Femoral D (Universal)	60	64
Femoral E (Universal)	62	66
Femoral F (Universal)	66	70
Femoral G (Universal)	70	74
Femoral H (Universal)	74	78

COMPONENT	AP	ML
Tibial Base 1 (Metalback)	40	59
Tibial Base 2 (Metalback)	40	62
Tibial Base 3 (Metalback)	42	66
Tibial Base 4 (Metalback)	46	66
Tibial Base 5 (Metalback)	48	71
Tibial Base 6 (Metalback)	50	72
Tibial Base 7 (Metalback)	52	76
Tibial Base 8 (Metalback)	54	78

TIBIAL TRIAL LINER	
Purple	A1-2
Yellow	B1-2
Rust	C1-2
Red	C3-4
Gray	D1-2
Green	D3-4
Light blue	E3-4
Blue	E5-6
Magenta	F3-4
Brown	F5-6
Bright Yellow	G5-6
Orange	G7-8
Black	H7-8

PATELLA TRIAL IMPLANT	
Yellow	25
Green	28
Light blue	31
Brown	34

Freedom Total Knee [®] System (CR and PS)		
	 Maxx Orthopedics, Inc. Plymouth Meeting, PA 19462 USA	
Carefully read all instructions and be familiar with the surgical techniques prior to use.		

Please see the package insert for complete device description, product selection information, indications, contraindications, precautions, adverse effects, warnings, materials, sterilization and patient guidance associated with the Freedom Total Knee[®] System.

CAUTION: THIS DEVICE IS RESTRICTED TO SALE BY OR ON THE ORDER OF A LICENSED PHYSICIAN

WARNING: THIS DEVICE IS INTENDED FOR CEMENTED USE ONLY

Maxx Orthopedics, Inc. is the manufacturer of the Freedom Total Knee[®] System.

Freedom[®] Knee is distributed by Maxx Health, Inc. Freedom[®] Knee is manufactured by Maxx Orthopedics, Inc. FREEDOM and FREEDOM TOTAL KNEE are Registered Trademarks of Maxx Orthopedics, Inc.

Maxx Orthopedics, Inc. Brochure No. MXO-MP00014-R01 © Maxx Orthopedics, Inc. All rights reserved.